

Anti-NCAM antibody [MEM-188] {Tri-Color} (STJ16100192)

STJ16100192

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD56 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-188
Concentration	
Conjugation	Tri-Color
Purification	
Dilution Range	
Formulation	PBS with 4mg/ml BSA, 10% sucrose and 0.1% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	