

## Anti-MyeloHRP antibody [8F4] (STJ16100195)

STJ16100195

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-MyeloHRP is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA and Immunofluorescence research applications.
<b>Applications</b>	IHC-F, FC, ELISA, IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	8F4
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	FC 1:50 IHC 1:50 ELISA 1:50 WB 1:50 IF 1:50
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen**  
**Region**  
**Specificity**  
**Immunogen**  
**Sequence**