

Anti-HLA-G antibody [MEM-G/9] {Biotin} (STJ16100197)

STJ16100197

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-G is suitable for use in Flow Cytometry, Immunoprecipitation, Immunohistochemistry, Immunocytochemistry and ELISA research applications.
Applications	FC, IP, IHC-F, ICC, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-G/9
Concentration	
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	3135
Gene Symbol	HLA-G
Uniprot ID	HLAG_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	