

Anti-HLA-E antibody [MEM-E/02] (STJ16100201)

STJ16100201

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-E is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB, IHC-P
Host/Source	Mouse
Reactivity	Human, Primates

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-E/02
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3133
Gene Symbol	HLA-E
Uniprot ID	HLAE_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	