

Anti-Csk binding protein antibody [PAG-C1] (STJ16100205)

STJ16100205

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Csk binding protein is suitable for use in Immunoprecipitation, Western Blot and Immunohistochemistry research applications.
Applications	IP, WB, IHC-P
Host/Source	Mouse
Reactivity	Rat, Cow, Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PAG-C1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	55824
Gene Symbol	PAG1
Uniprot ID	PHAG1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081