

Anti-HLA-G antibody [MEM-G/11] (STJ16100210)

STJ16100210

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-G is suitable for use in Immunoprecipitation, Flow Cytometry, Immunohistochemistry and ELISA research applications.
Applications	IP, FC, IHC-F, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MEM-G/11
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	3135
Gene Symbol	HLA-G
Uniprot ID	HLAG_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	