

## Anti-Nectin-2 antibody [502-57] (STJ16100216)

STJ16100216

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-Nectin-2 is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	IHC-F, FC, ELISA, IP, WB
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse, Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	502-57
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	FC 1:50 IHC 1:50 ELISA 1:50 WB 1:50 IP 1:50
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG2a
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">5819</a>
<b>Gene Symbol</b>	<a href="#">NECTIN2</a>
<b>Uniprot ID</b>	<a href="#">NECT2_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	