

Anti-Neurofilament antibody [RNF403] (STJ16100223)

STJ16100223

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Neurofilament is suitable for use in Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	IF, IHC-F, IHC-p, WB
Host/Source	Mouse
Reactivity	Human, Hamster, Monkey, Rat, Xenopus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RNF403
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB 1:10-1:20 IHC 1:10-1:20
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

