

Anti-Aflatoxin B1 & B2 antibody [AFT14/CBL03] (STJ16100224)

STJ16100224

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Aflatoxin B1 & B2 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Mouse
Reactivity	Aspergillus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AFT14/CBL03
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

