

## Anti-Basal cell Cytokeratin antibody [RCK103] (STJ16100229)

STJ16100229

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Basal cell Cytokeratin is suitable for use in Flow Cytometry, Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	FC, ICC, IHC-F, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Canine, Chicken, Guinea Pig, Hamster, Quail, Rabbit, Rat, Swine, Zebrafish

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	RCK103
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	FC (1:00-1:200) ICC IHC-F (1:100-1:200) WB (1:100-1:1000)
<b>Formulation</b>	PBS with 0.09% sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	

