

Anti-Keratin 18 antibody [M9] (STJ16100235)

STJ16100235

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 18 is suitable for use in Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	IHC-F, WB, FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M9
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F 1:150 IHC-P 1:10
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3875
Gene Symbol	KRT18
Uniprot ID	K1C18_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081