

Anti-Keratin 5+8 antibody [RCK102] (STJ16100236)

STJ16100236

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 5+8 is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-F, IHC-P, IF, WB
Host/Source	Mouse
Reactivity	Dog, Hamster, Human, Mouse, Pig, Rabbit, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RCK102
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:5-1:10
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence