

## Anti-Glial Fibrillary Acidic Protein antibody [GF-01] (STJ16100238)

STJ16100238

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Glial Fibrillary Acidic Protein is suitable for use in Immunohistochemistry, Immunoprecipitation, Western Blot and Immunocytochemistry research applications.
<b>Applications</b>	IHC-P, IP, WB, ICC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Pig, Cat, Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	GF-01
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">2670</a>
<b>Gene Symbol</b>	<a href="#">GFAP</a>
<b>Uniprot ID</b>	<a href="#">GFAP_HUMAN</a>
<b>Immunogen</b>	
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081