

Anti-Keratin 7 antibody [OVTL12/30] (STJ16100239)

STJ16100239

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 7 is suitable for use in research applications.
Applications	
Host/Source	Mouse
Reactivity	Universal

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OVTL12/30
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3855
Gene Symbol	KRT7
Uniprot ID	K2C7_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081