

Anti-Keratin 19 antibody [RCK108] (STJ16100245)

STJ16100245

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 19 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, IHC-p, WB
Host/Source	Mouse
Reactivity	Human, Zebrafish

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RCK108
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-F IHC-p (1:10-1:50) WB
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3880
Gene Symbol	KRT19
Uniprot ID	K1C19_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081