

Anti-Keratin 5 and 18 antibody [C-50] (STJ16100253)

STJ16100253

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 5 and 18 is suitable for use in Immunohistochemistry, Immunoprecipitation, Western Blot and Immunocytochemistry research applications.
Applications	IHC-P, IP, WB, ICC
Host/Source	Mouse
Reactivity	Rat, Hamster, Pig, Sheep, Dog, Cow, Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	C-50
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081