

Anti-Neurofilament antibody [403] (STJ16100255)

STJ16100255

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Neurofilament is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-F, IHC-P, IF, WB
Host/Source	Mouse
Reactivity	Bovine, Guinea Pig, Human, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	403
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ICC 1:10 IF 1:10 IHC 1:10
Formulation	Contains 0.1% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	