

## Anti-Epithelial Keratin Type II antibody [AE-3] (STJ16100262)

STJ16100262

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Epithelial Keratin Type II is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC-P, IHC-F, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rabbit, Rat, Bovine, Chicken

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	AE-3
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	WB 1:100 IHC 1:150-1:200
<b>Formulation</b>	PBS with 0.1% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	