

## Anti-Lewis B antibody [2-25LE] (STJ16100268)

STJ16100268

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Lewis B is suitable for use in ELISA, Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
<b>Applications</b>	ELISA, FC, IHC-F, IHC-P, IF, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human, Guinea Pig

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	2-25LE
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	FC (0.5-1.0 µg/million cells in 0.1 ml) WB (1.0-2.0 µg/ml) IF (0.5-1 0 µg/ml) ELISA (tracer: 0.001-100 µg/ml) IHC frozen and paraffin sections (1-2 µg/ml for 30 min at RT).
<b>Formulation</b>	PBS with 0.02% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen**  
**Region**  
**Specificity**  
**Immunogen**  
**Sequence**

