

Anti-Histone H1 antibody [AE-4] (STJ16100269)

STJ16100269

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Histone H1 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC, IHC-F, IHC-P, IF
Host/Source	Mouse
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AE-4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC (1.0 µg – 2.0 µg/million cells in 0.1 ml) IF (0 5-1 0 µg/ml) IHC frozen paraffin sections (1-2 µg/ml for 30 min at RT).
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

