

Anti-Lewis A antibody [7LE] (STJ16100271)

STJ16100271

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lewis A is suitable for use in ELISA, Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	ELISA, FC, IHC-F, IHC-P, IF, WB
Host/Source	Mouse
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	7LE
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC (0.5-1.0 µg/million cells in 0.1 ml) WB (1.0-2.0 µg/ml) IF (0.5-1 0 µg/ml) ELISA (tracer: 0.001-100 µg/ml) IHC frozen and paraffin sections (1-2 µg/ml for 30 min at RT).
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

