

Anti-JAM-1 antibody [M.Ab.F11] (STJ16100273)

STJ16100273

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-JAM-1 is suitable for use in Flow Cytometry, ELISA, Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	FC, ELISA, IF, IP, IHC-P, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M.Ab.F11
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:50 IP 1:50 IHC 1:50 WB 1:50 IP 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	50848
Gene Symbol	F11R
Uniprot ID	JAM1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	