

Anti-Desmoglein-2 antibody [6D8] (STJ16100276)

STJ16100276

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Desmoglein-2 is suitable for use in Immunohistochemistry, Immunofluorescence, Immunoprecipitation and Western Blot research applications.
Applications	IHC-F, IF, IP, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6D8
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:10 WB 1:10 IP 1:10 IF 1:10 IHC 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	1829
Gene Symbol	DSG2
Uniprot ID	DSG2_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	