

Anti-Desmoglein-3 antibody [5G11] (STJ16100277)

STJ16100277

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Desmoglein-3 is suitable for use in Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	IF, IP, IHC-P, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5G11
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	WB 1:10 IP 1:10 IF 1:10 IHC 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	1830
Gene Symbol	DSG3
Uniprot ID	DSG3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	