

Anti-Calponin-1 antibody [EP798Y] (STJ16100318)

STJ16100318

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Calponin-1 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EP798Y
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:50-1:200
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	1264
Gene Symbol	CNN1
Uniprot ID	CNN1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

