

Anti-CD14 antibody [EPR3653] (STJ16100322)

STJ16100322

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-CD14 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EPR3653
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:50-1:200
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	929
Gene Symbol	CD14
Uniprot ID	CD14_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

