

Anti-CD2 antibody [MRQ-11] (STJ16100326)

STJ16100326

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD2 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MRQ-11
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:10-1:50
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	914
Gene Symbol	CD2
Uniprot ID	CD2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

