

Anti-Cytokeratin 14 antibody [LL002] (STJ16100349)

STJ16100349

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cytokeratin 14 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LL002
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1: 100 – 1:500
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG3
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3861
Gene Symbol	KRT14
Uniprot ID	K1C14_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

