

Anti-Cytokeratin 5 antibody [EP1601Y] (STJ16100352)

STJ16100352

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Cytokeratin 5 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EP1601Y
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1: 50 – 1:200
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3852
Gene Symbol	KRT5
Uniprot ID	K2C5_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

