

Anti-Cytokeratin antibody [AE1/AE3] (STJ16100357)

STJ16100357

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Cytokeratin is suitable for use in Immunohistochemistry research applications. |
| Applications | IHC-P |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | AE1/AE3 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | IHC 1: 100 – 1:500 |
| Formulation | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Isotype | IgG1k |
| Storage Instruction | Store at 2-8°C upon receipt. |

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence

