

## Anti-Cytokeratin antibody [OSCAR] (STJ16100358)

STJ16100358

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Cytokeratin is suitable for use in Immunohistochemistry research applications. |
| <b>Applications</b>      | IHC-P   |
| <b>Host/Source</b>       | Mouse   |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Monoclonal  |
| <b>Clone ID</b>            | OSCAR   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      | IHC 1: 25 – 1:100   |
| <b>Formulation</b>         | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| <b>Isotype</b>             | IgG2a   |
| <b>Storage Instruction</b> | Store at 2-8°C upon receipt.                                      |

### TARGET INFORMATION

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen Region**  
**Specificity**  
**Immunogen Sequence**

