

Anti-Ep-CAM/Epithelial Specific Antigen antibody [Ber-EP4] (STJ16100363)

STJ16100363

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Ep-CAM/Epithelial Specific antigen is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Ber-EP4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:50 – 1:200
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4072
Gene Symbol	EPCAM
Uniprot ID	EPCAM_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	