

Anti-Glycophorin A antibody [GA-R2 (HIR2)] (STJ16100373)

STJ16100373

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Glycophorin A is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	GA-R2 (HIR2)
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:100 – 1:500
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG2bk
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2993
Gene Symbol	GYPA
Uniprot ID	GLPA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

