

Anti-p21 antibody [DCS-60.2] (STJ16100407)

STJ16100407

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-p21 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DCS-60.2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:50 – 1:200
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	1026
Gene Symbol	CDKN1A
Uniprot ID	CDN1A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

