

Anti-Thrombomodulin antibody [1009] (STJ16100458)

STJ16100458

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Thrombomodulin is suitable for use in Immunohistochemistry research applications.
Applications	IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1009
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:25 – 1:100
Formulation	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7056
Gene Symbol	THBD
Uniprot ID	TRBM_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

