

Anti-Myosin antibody [414] (STJ16100469)

STJ16100469

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Myosin is suitable for use in Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	IHC-F, IHC-P, IF, WB
Host/Source	Mouse
Reactivity	Chicken, Human, Mouse, Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	414
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IF 1:10 IHC 1:5-1:10
Formulation	Contains 0.09% sodium azide
Isotype	IgM
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	