

Anti-Collagen IV antibody [MC4-HA] (STJ16100470)

STJ16100470

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Collagen IV is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MC4-HA
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ELISA 1:1000 IHC 1:20
Formulation	Purified and lyophilized
Isotype	IgM
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	