

Anti-Beta-Tubulin antibody [TU-06] (STJ16100474)

STJ16100474

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Beta-Tubulin is suitable for use in Western Blot, Immunohistochemistry and Immunocytochemistry research applications.
Applications	WB, IHC-P, ICC
Host/Source	Mouse
Reactivity	A.thaliana, Mouse, Rat, Pig, Nicotiana, Paramecium, Chicken, Human, Fish, Tetrahymena

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TU-06
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgM
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081