

Anti-Alpha-Tubulin antibody [TU-02] (STJ16100478)

STJ16100478

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Alpha-Tubulin is suitable for use in ELISA and Immunohistochemistry research applications.
Applications	ELISA, IHC-F
Host/Source	Mouse
Reactivity	Human, Mouse, Rat, Pig, Plants

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TU-02
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ELISA (solid phase: not known; tracer: 0.001-100 µg/ml) IHC paraffin sections (1-2 µg/ml for 30 min at RT).
Formulation	PBS with 0.02% sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	

