

Anti-JAM-C antibody [CRAM-18 F26] (STJ16100484)

STJ16100484

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-JAM-C is suitable for use in Immunohistochemistry, Flow Cytometry, ELISA and Immunoprecipitation research applications.
Applications	IHC-F, FC, ELISA, IP
Host/Source	Mouse
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CRAM-18 F26
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	83700
Gene Symbol	JAM3
Uniprot ID	JAM3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	