

Anti-TLR2 antibody [mT2.7] {Biotin} (STJ16100486)

STJ16100486

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR2 is suitable for use in research applications.
Applications	
Host/Source	Mouse
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	mT2.7
Concentration	
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7097
Gene Symbol	TLR2
Uniprot ID	TLR2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	