

Anti-CD41b antibody [HIP2] (STJ16100488)

STJ16100488

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD41b is suitable for use in Flow Cytometry and Immunohistochemistry research applications.
Applications	FC, IHC-F
Host/Source	Mouse
Reactivity	Primates, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HIP2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgG3
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	