

Anti-TLR2 antibody [TL2.1] (STJ16100490)

STJ16100490

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TLR2 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, FC, IF, IP, IHC-P, WB
Host/Source	Mouse
Reactivity	Human, Canine, Cynomolgus, Monkey, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	TL2.1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7097
Gene Symbol	TLR2
Uniprot ID	TLR2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	