

Anti-TLR4/MD2 antibody [MTS510] (STJ16100494)

STJ16100494

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-TLR4/MD2 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunofluorescence and Immunoprecipitation research applications.
Applications	IHC-F, FC, IF, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MTS510
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	23643
Gene Symbol	LY96
Uniprot ID	LY96_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	