

Anti-TLR4/MD2 antibody [MTS510] {Biotin} (STJ16100495)

STJ16100495

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-TLR4/MD2 is suitable for use in research applications.
Applications	
Host/Source	Rat
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MTS510
Concentration	
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	23643
Gene Symbol	LY96
Uniprot ID	LY96_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	