

Anti-CD32 antibody [AT10] (STJ16100502)

STJ16100502

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD32 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	FC, IHC-F, IP
Host/Source	Mouse
Reactivity	Human, Pig, Rhesus Monkey, Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AT10
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	TBS
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	2212
Gene Symbol	FCGR2A
Uniprot ID	FCG2A_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081