

Anti-HIV p24 antibody [ND1] (STJ16100518)

STJ16100518

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HIV p24 is suitable for use in Western Blot, Immunocytochemistry and ELISA research applications.
Applications	WB, ICC, ELISA
Host/Source	Mouse
Reactivity	Human, Immunodeficiency, Virus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ND1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	