

Anti-CD89 antibody [MIP8a] {PE} (STJ16100534)

STJ16100534

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD89 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MIP8a
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2204
Gene Symbol	FCAR
Uniprot ID	FCAR_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	