

Anti-CD91 antibody [A2Mr alpha-2] (STJ16100538)

STJ16100538

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD91 is suitable for use in Western Blot, Immunoprecipitation, Flow Cytometry, Immunohistochemistry and ELISA research applications.
Applications	WB, IP, FC, IHC-F, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A2Mr alpha-2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence