

Anti-CD91 antibody [A2Mr alpha-2] {PE} (STJ16100539)

STJ16100539

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD91 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A2Mr alpha-2
Concentration	
Conjugation	PE
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide and 1% BSA, 5% sucrose
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	